# lumie® Bodyclock™

STARTER 30

**OPERATING INSTRUCTIONS**

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**BRUKSANVISNING**

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# OPERATING INSTRUCTIONS

## Bodyclock Starter 30

---

- 1 Safety
- 2 Getting started
- 3 Programming
- 4 User guide
- 5 Problem-solving

## INTRODUCTION

Bodyclock Starter wakes you with a sunrise that naturally helps to regulate your sleep/wake pattern. You'll wake in a better mood, find it easier to get up and feel more energetic during the day. In clinical trials, 'dawn simulators' like Bodyclock Starter have been shown to improve mood, energy, productivity and the quality of sleep and awakening\*. It can also help SAD (Seasonal Affective Disorder) and winter blues sufferers cope with dark winter mornings.

Please read these instructions carefully to get the most out of Bodyclock Starter.

\*Acta Psychiatr Scand 1993 Jul; 88 (1):67-71


## 1 Safety

- Keep the lamp cover over the bulb. Remove only when cold and when the lamp is unplugged.
- The lamp cover will be warm to the touch. Avoid covering it and keep clear anything that may be damaged by heat.
- Don't open the casing for any reason. There are no user-serviceable parts inside.
- Keep away from water and damp. For indoor use only.
- Maximum 42W halogen bulb/40W candle bulb.

## 2 Getting started

The best position for Bodyclock Starter is by your bedside.

- Remove the lamp cover with a gentle upward tug.
- Screw in the bulb firmly and replace the cover.
- Plug into the mains and you're ready for *Programming*.

### 3 Programming

Just a few programs provide you with the basic Bodyclock functions: a 30-minute sunrise to wake you up; a back-up alarm beep for added reassurance; an optional 30-minute sunset to help you sleep.

**Clock time** (default 13:00) – 5EŁc shows on the display

**Alarm time** (default 07:00 – 5EŁR shows on the display/\* flashes

**12/24hr clock** (default 24hr) – 12hr or 24hr shows on the display/  
**am pm** flashes

**Audible alarm** (default OFF) – 0FF or 0n shows on the display/🔊 flashes

- Press and hold Ⓞ button until e.g. 5EŁc shows on display.
- Use –/+ buttons to change settings.
- Press Ⓞ button to skip to the next program.
- Wait 12 seconds for the display to revert to normal or quick press on Ⓞ to skip to next program.

All settings are saved for around an hour if the unit is unplugged or if there's a power failure.

## 4 User guide

### **Sunset and sunrise – recommended routine**

Use Bodyclock Starter to wake up and go to sleep at the same time every day.

- Use the -/+ buttons to turn on, adjust and turn off the light.
- When you're ready to sleep, press  to turn on the sunrise alarm (alarm time shows/ displayed). Display auto-dim is triggered when the light is at low level and the alarm is set.
- This automatically starts the sunset and the light will slowly fade down.
- If you *don't* want a sunset, simply turn off the light when you're ready.
- When you wake up, press  to turn off the alarm.

*Note: Both the audible alarm and the light will turn off after 20 minutes if they are not turned off manually.*

## 5 Problem-solving

Replacement bulbs are available from Lumie.

- **I ran a test sunrise but the light came up really quickly.**

When setting a test alarm make sure you allow enough time for a full sunrise – the default sunrise takes 30 minutes. For example, if the actual time is 11:00 the alarm time must be set no earlier than 11:30, otherwise the sunrise speeds up in order to finish on time.

- **The light will not come on but the clock still works.**

Check to make sure the bulb is fully screwed into the socket. If this doesn't work, try replacing the bulb.

- **Neither the clock nor the light works, or the clock is not doing what I expect it to.**

Turn off the power at the wall. Press and hold  as you turn the power back on. This re-sets all the program options to their default state. If this doesn't work, check the fuse in the plug.

Please call the Lumie Careline (01954 780500) if you have any other problems.

## Warranty

This unit is guaranteed against any manufacturing defect for 24 months from the date of purchase. This is in addition to your statutory rights. If you are returning your Bodyclock Starter for service, it must be in its original packaging including internal pieces. Please call us for further advice.

## Technical specification

Specification subject to change without notice.

- Lumie Bodyclock Starter 30
- 230V 50Hz 50W
- Lamp E14 (SES) 35mm halogen capsule 240V 42W MAX
- Operating temperature 0°C – 35°C
- Class I Medical Device

## Contact us

Lumie is a trademark of Outside In, Europe's leading specialist in light therapy. Our products are based on the latest scientific studies and we focus solely on light therapy and its benefits. If you have any suggestions that you feel could help us improve our products or service, or if there are problems you think we should know about, please call the Lumie Careline on 01954 780500 or email [info@lumie.com](mailto:info@lumie.com).

Lumie

3 The Links, Trafalgar Way

Bar Hill

Cambridge

CB23 8UD

UK

Lumie is a trademark of Outside In (Cambridge) Limited  
incorporated in England and Wales.

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© Lumie 2009
Lumie is een handelsmerk van Outside In (Cambridge) Limited